

Force Tensiometer

K6

Specifications



Product group specifications

K6

Force measurement

Calibration	external weight
Calibration weight	CP11 or CP12 (optional)

Sample stage

Travel distance	85 mm
Simple platform	yes

Measurement specifications

K6

Du Noüy ring

Results	surface tension (SFT)/interfacial tension (IFT)
Range	1 to 90 mN/m
Resolution	0.5 mN/m
Correction method	linear correction

General specifications

K6

Housing and peripherals

Built-in bubble level	yes
-----------------------	-----

Instrument dimensions

Footprint	230 mm × 270 mm (W × D)
Height	315 mm
Weight (without accessories)	3 kg